

HUMAN TYRO3 PROTEIN, HIS TAG

货号: 11270

产品全名: 人 TYRO3 蛋白

规格: 10/50/100 µg

基因符号 BYK;Dtk;Etk-2;Rek;RSE;Sky;Tif

目标蛋白: TYRO3

UNIPROT ID: Q06418

描述: Recombinant human TYRO3 protein with C-terminal 6xHis tag

背景: The gene is part of a 3-member transmembrane receptor kinase receptor family with a processed pseudogene distal on chromosome 15. The encoded protein is activated by the products of the growth arrest-specific gene 6 and protein S genes and is involved in controlling cell survival and proliferation, spermatogenesis, immunoregulation and phagocytosis. The encoded protein has also been identified as a cell entry factor for Ebola and Marburg viruses. [provided by RefSeq, May 2010]

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 42.4 kDa after removal of the signal peptide. The apparent molecular mass of TYRO3-His is approximately 55-100 kDa due to glycosylation.

分子特征: TYRO3(Ala41-Ser428) 6×His tag

纯化: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Human TYRO3 Protein, His Tag on SDS-PAGE under reducing condition.