

## **IGF1R (DMC684) IGG1 CHIMERIC MAB**

货号: 28282

产品全名: IGF1R(DMC684) IgG1 Chimeric 单抗

基因符号 CD221; IGF1R; IGF1; JTK13

描述: IGF1R antibody(DMC684) IgG1 Chimeric 单抗

背景: This receptor binds insulin-like growth factor with a high affinity. It has tyrosine kinase activity. The insulin-like growth factor I receptor plays a critical role in transformation events. Cleavage of the precursor generates alpha and beta subunits. It is highly overexpressed in most malignant tissues where it functions as an anti-apoptotic agent by enhancing cell survival. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq; May 2014]

经过测试的应用: Flow Cyt

推荐稀释比: Flow Cyt 1:100

种属反应性: Rabbit

亚型: Rabbit:Human Fc chimeric IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 IGF1R

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. IGF1R protein is highly expressed on the surface of Expi293 cell membrane. Flow cytometry analysis with Anti-IGF1R (DMC684) on Expi293 cells transfected with human IGF1R (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram), and Isotype antibody on Expi293 transfected with irrelevant protein (Orange histogram).