

## IL7RA (DMC272) IGG1 CHIMERIC MAB

货号: 28180

产品全名: IL7RA(DMC272) IgG1 Chimeric 单抗

基因符号 CD127; CDW127; IL-7R-alpha; IL7RA; ILRA

描述: IL7RA antibody(DMC272) IgG1 Chimeric 单抗

背景: The protein encoded by this gene is a receptor for interleukin 7 (IL7). The function of this receptor requires the interleukin 2 receptor; gamma chain (IL2RG); which is a common gamma chain shared by the receptors of various cytokines; including interleukins 2; 4; 7; 9; and 15. This protein has been shown to play a critical role in V(D)J recombination during lymphocyte development. Defects in this gene may be associated with severe combined immunodeficiency (SCID). Alternatively spliced transcript variants have been found.

经过测试的应用: Flow Cyt

推荐稀释比: Flow Cyt 1:100

种属反应性: Rabbit

亚型: Rabbit:Human Fc chimeric IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 IL7RA

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with Anti-IL7RA (DMC272) on Expi293 cells transfected with human IL7RA (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).