

CD138 MAB (INDATUXIMAB)

货号: 28023

产品全名: CD138 (Indatuximab) 单抗

基因符号 SDC1;Syndecan-1;CD138;SYND1;SDC

描述: CD138 (indatuximab) 单抗

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释比: Flow Cyt 1:100

种属反应性: Chimeric

亚型: IgG4

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 CD138

成分: Lyophilized from sterile PBS, pH 7.4. 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. HEK293 cell line transfected with irrelevant protein (B) and human CD138 (A) were surface stained with anti-CD138 neutralizing antibody $1\mu\text{g}/\text{ml}$ (indatuximab) followed by Alexa 488-conjugated anti-human IgG secondary antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
