

**HUMAN KCNK9 FULL LENGTH PROTEIN**

货号: 11044

产品全名: 人 KCNK9 Full Length 蛋白

规格: 10/50/100 µg

基因符号 BIBARS; K2p9.1; KT3.2; TASK-3; TASK3; TASK32

目标蛋白: KCNK9

**UNIPROT ID:** Q9NPC2

描述: 人 KCNK9 full length protein-synthetic nanodisc

背景: This gene encodes a protein that contains multiple transmembrane regions and two pore-forming P domains and functions as a pH-dependent potassium channel. Amplification and overexpression of this gene have been observed in several types of human carcinomas. This gene is imprinted in the brain, with preferential expression from the maternal allele. A mutation in this gene was associated with Birk-Barel dysmorphism syndrome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2017]

物种/宿主: HEK293

分子量: The human full length KCNK9 protein has a MW of 42.3 kDa

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.