

## MOUSE BMP6 PROTEIN, HIS TAG

货号: 12197

产品全名: Mouse BMP6 蛋白

规格: 10/50/100 µg

基因符号 BMP-6;VGR-1

目标蛋白: BMP6

**UNIPROT ID:** P20722

**描述:** Recombinant mouse BMP6 protein with C-terminal 6xHis tag

**背景:** This gene encodes a secreted ligand of the TGF-beta (transforming growth factor-beta) superfamily of proteins. Ligands of this family bind various TGF-beta receptors leading to recruitment and activation of SMAD family transcription factors that regulate gene expression. The encoded preproprotein is proteolytically processed to generate each subunit of the disulfide-linked homodimer. This protein regulates a wide range of biological processes including iron homeostasis, fat and bone development, and ovulation. Mice lacking this gene exhibit delayed ossification of the sternum, iron overload, and reduced fertility in females. [provided by RefSeq, Jul 2016]

**物种/宿主:** HEK293

**分子量:** The protein has a predicted molecular mass of 54.9 kDa after removal of the signal peptide. The apparent molecular mass of mBMP6-His is approximately 35-70 kDa due to glycosylation.

**分子特征:** Mouse BMP6(Cys21-His510) 6×His tag

**纯化:** The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Figure 1. Mouse BMP6 Protein, His Tag on SDS-PAGE under reducing condition.