

MOUSE CCR8 PROTEIN, HFC TAG

货号: 12158

产品全名: Mouse CCR8 蛋白

规格: 10/50/100 µg

基因符号 CC-

CKR-8;CCR-8;CDw198;CKRL1;CMKBR8;CMKBRL2;CY6;GPRCY6;TER1

目标蛋白: CCR8

UNIPROT ID: P56484

描述: Recombinant mouse CCR8 protein with C-terminal human Fc tag

背景: Receptor for the CCL1/SCY1/TCA-3 chemokine.[UniProtKB/Swiss-Prot Function]

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 30.0 kDa after removal of the signal peptide. The apparent molecular mass of mCCR8-hFc is approximately 35-55 kDa due to glycosylation.

分子特征: Mouse CCR8(Met1-Met33) hFc(Glu99-Ala330)

纯化: The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Mouse CCR8 Protein, hFc Tag on SDS-PAGE under reducing condition.