

## MOUSE CD138 PROTEIN, HFC TAG

货号: 12138

产品全名: Mouse CD138 蛋白

规格: 10/50/100 µg

基因符号 AA408134;AA409076;CD138;Sstn;syn-1;Synd;Synd1

目标蛋白: CD138

**UNIPROT ID:** P18828

**描述:** Recombinant mouse CD138 protein with C-terminal human Fc tag

**背景:** Cell surface proteoglycan that bears both heparan sulfate and chondroitin sulfate and that links the cytoskeleton to the interstitial matrix. Regulates exosome biogenesis in concert with SDCBP and PDCD6IP.[UniProtKB/Swiss-Prot Function]

**物种/宿主:** HEK293

**分子量:** The protein has a predicted molecular mass of 50.3 kDa after removal of the signal peptide. The apparent molecular mass of mCD138-hFc is approximately 55-70 kDa due to glycosylation.

**分子特征:** Mouse CD138(Gln23-Leu254) hFc(Glu99-Ala330)

**纯化:** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Mouse CD138 Protein, hFc Tag on SDS-PAGE under reducing condition.