

## MOUSE CD27 PROTEIN, HFC TAG

货号: 12135

产品全名: Mouse CD27 蛋白

规格: 10/50/100 µg

基因符号 S152;Tnfrsf7;Tp55

目标蛋白: CD27

**UNIPROT ID:** P41272

**描述:** Recombinant mouse CD27 protein with C-terminal human Fc tag

**背景:** Receptor for CD70/CD27L. May play a role in survival of activated T-cells. May play a role in apoptosis through association with SIVA1.[UniProtKB/Swiss-Prot Function]

**物种/宿主:** HEK293

**分子量:** The protein has a predicted molecular mass of 43.9 kDa after removal of the signal peptide. The apparent molecular mass of mCD27-hFc is approximately 55-70 kDa due to glycosylation.

**分子特征:** Mouse CD27(Pro24-Arg182) hFc(Glu99-Ala330)

**纯化:** The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization.

**储存和运输:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.



Figure1. Mouse CD27 Protein, hFc Tag on SDS-PAGE under reducing condition.