

MOUSE MESOTHELIN PROTEIN, HIS TAG

货号: 12164

产品全名: Mouse Mesothelin 蛋白

规格: 10/50/100 µg

基因符号 Mesothelin;Pre-pro-megakaryocyte-potentiating factor

目标蛋白: Mesothelin

UNIPROT ID: Q61468

描述: Recombinant mouse Mesothelin protein with C-terminal 6xHis tag

背景: This gene encodes a preproprotein that is proteolytically processed to generate two protein products, megakaryocyte potentiating factor and mesothelin. Megakaryocyte potentiating factor functions as a cytokine that can stimulate colony formation of bone marrow megakaryocytes. Mesothelin is a glycosylphosphatidylinositol-anchored cell-surface protein that may function as a cell adhesion protein. This protein is overexpressed in epithelial mesotheliomas, ovarian cancers and in specific squamous cell carcinomas. Alternative splicing results in multiple transcript variants, at least one of which encodes an isoform that is proteolytically processed. [provided by RefSeq, Feb 2016]

物种/宿主: HEK293

分子量: The protein has a predicted molecular mass of 35.0 kDa after removal of the signal peptide. The apparent molecular mass of mMesothelin-His is approximately 35-40 kDa due to glycosylation.

分子特征: Mouse Mesothelin(Asp298-Ser600) 6×His tag

纯化: The purity of the protein is greater than 85% as determined by SDS-PAGE and Coomassie blue staining.

Formulation & Reconstitution: Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.




Figure 1. Mouse Mesothelin Protein, His Tag on SDS-PAGE under reducing condition.