

MUC1 PAB

产品名称: MUC1 兔多抗

货号: 21123

基因符号: MUC1

描述: MUC1 兔多抗

背景: Mucin-1 proteins line the apical surface of highly polarized epithelial cells in the lungs, stomach, intestines, and eyes. Mucins protect the body from infection by pathogen binding to oligosaccharides in the extracellular domain. Overexpression of MUC1 is often associated with colon, breast, ovarian, lung, and pancreatic cancers.

免疫原: A synthetic peptide from the internal region of MUC1 of human origin.

经过测试的应用: ELISA, WB, IHC, IP

推荐稀释比:

ELISA: 1:100-1:5000

WB: 1:100-1:1000

IHC: 1:50-1:100

Concentration: 0.2 mg/ml

种属反应性: Rabbit

形式: Liquid

克隆性: Polyclonal


纯化: Purified from serum

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: MUC1 antibody recognizes MUC1 of vertebrates.

储存条件: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

蛋白质印迹:



Whole-cell lysates of MEF (lane 1), HEK293 (lane 2) and NIH3T3 (lane 3) cells were used for western blot with the polyclonal Anti-MUC1 antibody (Cat # 21123).