

## RAB35-GTP

产品名称: Rab35-GTP

货号: 26922

规格: 30  $\mu$ L

基因符号: RAB35

描述: Rab35-GTP 单抗

背景: The Rab family of GTPases regulates eukaryotic vesicular membrane traffic. Rab35 is with both plasma membrane and endosomal localization, and has been implicated in diverse processes that include T-cell receptor recycling, oocyte yolk protein recycling and cytokinesis. Rab35 regulates neurite outgrowth in neuronal-like cells, and can induce protrusions.

免疫原: Recombinant full length protein of active Rab35

经过测试的应用: IP, IHC and IF (**Not applicable for WB since WB denatures the GTPase**)

推荐稀释比:

IP: 1  $\mu$ g for 1~2 mg total cellular proteins

IHC, IF: 1:50~1:250

**Concentration:** 1 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG

纯化: Purified from ascites


**Preservative:** No

成分: PBS (without  $Mg^{2+}$  and  $Ca^{2+}$ ), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Rab35-GTP monoclonal antibody recognizes active Rab35 proteins of vertebrates.

储存条件: Store at  $-20^{\circ}\text{C}$ . Avoid repeated freezing and thawing

免疫沉淀/蛋白质印迹:



IP: anti-active Rab35 mAb

IB: anti-Rab35 pAb

**Rab35 Activation Assay.** Purified full-length Rab35 proteins were immunoprecipitated after treated with GDP (lane 1) or GTP $\gamma$ S (lane 2). Immunoprecipitation was done with the anti-Rab35-GTP monoclonal antibody (货号: 26922). Immunoblot was with an anti-Rab35 兔多抗 (货号: 21078).