

## RAB3A MAB

产品名称: RAB3A 小鼠单抗

货号: 26031

基因符号: RAB3A

描述: Rab3A 小鼠单抗

背景: Rab3A is a member of the Rab-family GTPases and is specifically targeted to synaptic vesicles.

免疫原: Recombinant full length protein of Rab3A.

经过测试的应用: ELISA, WB, IHC, IP

推荐稀释比:

ELISA: 1:100-1:5000

WB: 1:100-1:2000

IHC: 1:100-1:1000

**Concentration:** 1 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG


纯化: Purified from ascites

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Rab3A antibody recognizes Rab3A of vertebrates.

储存条件: Store at -20°C. Avoid repeated freezing and thawing

免疫沉淀/蛋白质印迹:



**Wb:anti-Rab3A mAb**

Western blot analysis of recombinant GST-Rab3A proteins.