

RAB5 PAB

产品名称: Rab5 兔多抗

货号: 21151

基因符号: RAB5

描述: Rab5 兔多抗

背景: Rab5, consists of three isoforms, Rab5A, Rab5B, and Rab5C, is a small GTPase that is localized to early endosomes. It regulates the fusion between endocytic vesicles and early endosomes, as well as the homotypic fusion between early endosomes.

免疫原: Three synthetic peptides from the internal region of human Rab5A, Rab5B and Rab5C proteins, respectively

经过测试的应用: ELISA, WB, IHC, IP

推荐稀释比:

ELISA: 1:100-1:5000

WB: 1:100-1:2000

IHC: 1:100-1:1000

Concentration: 0.2 mg/ml

种属反应性: Rabbit

形式: Liquid

克隆性: Polyclonal

亚型: IgG


纯化: Purified from serum

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Rab5 antibody recognizes Rab5A, Rab5B and Rab5C from vertebrates.

储存条件: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

蛋白质印迹:



WB: anti-Rab5 pAb

EC-1 whole-cell lysate was used for Western blot with anti Rab9A 兔多抗 (cat # 21073).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
