

RAB7 MAB

产品名称: Rab7 小鼠单抗

货号: 26055

基因符号: RAB7

描述: Rab7 小鼠单抗

背景: Rab7 belongs to the Rab family of small GTPases. Rab7 is located in the late endosome, lysosome, cytoplasmic vesicle, phagosome and melanosome. It is involved in the late endocytic transport and contributes to the maturation of phagosomes.

免疫原: purified His-Rab7 蛋白

经过测试的应用: ELISA, WB, IHC, IP

推荐稀释比:

ELISA: 1:100–1:5000

WB: 1:100–1:2000

IHC: 1:100–1:1000

Concentration: 0.8 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG

纯化: Purified from ascites

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 59% glycerol

种属反应性: Rab7 antibody recognizes Rab7 proteins of vertebrates.

储存条件: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

蛋白质印迹:



Western blot analysis of recombinant Rab7 proteins. Purified His-tagged Rab7 protein (货号: 10167) was blotted with anti-Rab7 monoclonal antibody (Cat # 26055).