

SAR1 MAB

产品名称: Sar1 小鼠单抗

货号: 26021

基因符号: SAR1

描述: Sar1 小鼠单抗

背景: Small GTPases are a super-family of cellular signaling regulators. Sar1 is a member of these small GTPases. Sar1 functions as a molecular switch to control protein-protein and protein-lipid interactions that direct vesicle budding from the ER.

免疫原: Recombinant full length Sar1a protein.

经过测试的应用: ELISA, WB, IHC

推荐稀释比: ELISA: 1:100-1:5000WB: 1:100-1:2000IHC: 1:50-1:100

Concentration: 1 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG


纯化: Purified from ascites

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Sar1 antibody recognizes Sar1 proteins of vertebrates.

储存条件: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

Western blot:



蛋白质印迹: WB: anti-Sar1 mAb Western blot analysis of Sar1 proteins. EC-1

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC APPLICATIONS

whole-cell lysate was used for Western blots with anti-Sar1 小鼠单抗 (Cat# 26021).