

## IL17A MAB (SECUKINUMAB)

货号: 28131

产品全名: IL17A(Secukinumab) 单抗

基因符号 CTLA-8;TLA8;L-17;L-17A;L17

描述: IL17A(secukinumab) 单抗

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释比: Flow Cyt 1:100

种属反应性: Homo sapiens

亚型: IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 IL17A

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis under cell membrane permeable condition with 1  $\mu\text{g}/\text{mL}$  Anti-IL17A (secukinumab biosimilar) mAb (28131) on Expi293 cells transfected with Human IL17A protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).