

SEPTIN9 PAB

产品名称: Septin9 兔多抗

货号: 21126

基因符号: SEPTIN9

描述: Septin9 兔多抗

背景: Septin9 belongs to the septin family. It is essential for pseudopod protrusion, tumor cell migration and invasion.

免疫原: A synthetic peptide from the internal region of Septin9 isoform of human origin.

经过测试的应用: ELISA, WB, IHC, IP

推荐稀释比:

ELISA: 1:100-1:5000

WB: 1:100-1:1000

IHC: 1:50-1:100

Concentration: 0.2 mg/ml

种属反应性: Rabbit

形式: Liquid

克隆性: Polyclonal


纯化: Purified from serum

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Septin-9 antibody recognizes Septin-9 of vertebrates.

储存条件: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

蛋白质印迹:



WB: anti-Septin 9 pAb

Whole-cell lysate of HEK293T cells was used for western blot with the polyclonal Anti-Septin9 antibody (Cat # 21126).