

IL6 MAB (SILTUXIMAB)

货号: 28016

产品全名: IL6 (Siltuximab) 单抗

基因符号 IL6;Interleukin-6;BSF2;HSF;IFNB2

描述: IL6 (siltuximab) 单抗

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: ELISA; Flow Cyt

推荐稀释比: ELISA 1:5000-10000; Flow Cyt 1:100

种属反应性: Chimeric

亚型: IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 IL6

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form.

After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

Anti-IL6 (siltuximab biosimilar) mAb ELISA

0.2 μ g of Human IL6, mFc-His Tagged protein per well




Figure 1. ELISA plate pre-coated by 2 μ g/ml (100 μ l/well) Human IL6, mFc-His tagged protein 11160 can bind Anti-IL6 Neutralizing antibody in a linear range of 7.81-386.8 ng/ml.