

**TP53(R175H) PAB**

产品名称: TP53(R175H) 兔多抗

货号: 21161

基因符号: TP53(R175H)

描述: TP53(R175H) 兔多抗

背景: p53 protein acts as a tumor suppressor in many tumor types and induces growth arrest or apoptosis depending on the physiological circumstances and cell type. Defects in TP53 are found in many types of cancers.

免疫原: A synthetic peptide from the internal region of p53 which includes the mutation of R175H, human origin.

经过测试的应用: ELISA, WB, IHC, IP

推荐稀释比: ELISA: 1:100-1:5000WB: 1:100-1:2000IHC: 1:50-1:100

**Concentration:** 0.1 mg/ml

种属反应性: Rabbit

形式: : Liquid

克隆性: Polyclonal


亚型:

纯化: Purified from serum

成分: PBS (without  $Mg^{2+}$  and  $Ca^{2+}$ ), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: TP53(R175H) antibody recognizes TP53(R175H) of vertebrates.

储存条件: Store at  $-20^{\circ}C$ . Avoid repeated freezing and thawing

**蛋白质印迹:**

WB: Anti-TP53(R175H) pAb

Western blot analysis of recombinant TP53(R175H) and wild type p53 proteins.

Purified His-tagged TP53(R175H) protein (货号:10207) and wild type p53 protein (货

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC APPLICATIONS

号:10202) were blotted with Anti-TP53(R175H) rabbit antibody (货号: 21161).