

TP53(R282W) PAB

产品名称: TP53(R282W) 兔多抗

货号: 21165

基因符号: TP53

描述: TP53(R282W) 兔多抗

背景: p53 protein acts as a tumor suppressor in many tumor types and induces growth arrest or apoptosis depending on the physiological circumstances and cell type. Defects in TP53 are found in many types of cancers.

免疫原: A synthetic peptide near the C-terminus of human TP53 which includes the mutation of R282W.

经过测试的应用: ELISA, WB, IHC, IF

推荐稀释比: ELISA: 1:1000-1:5000 WB: 1:1000-1:5000 IHC: 1:100-1:1000

Concentration: 0.1 mg/ml

种属反应性: Rabbit

形式: Liquid

克隆性: Polyclonal


亚型:

纯化: Purified from serum

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: TP53(R282W) antibody recognizes TP53(R282W) of vertebrates.

储存条件: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing.



蛋白质印迹: WB: Anti-TP53(R282W) pAb Western blot analysis of recombinant TP53(R282W) and wild type proteins. Purified His-tagged TP53(R282W) protein (货号: 10204) and wild type protein (货号: 10202) were blotted with Anti-TP53(R282W) rabbit antibody (货号: 21165).

