

TRANSFERRIN MAB

产品名称: Transferrin Antibody

货号: 26050

基因符号: TF

描述: Transferrin 小鼠单抗

背景: Transferrins are iron-binding blood plasma glycoproteins that control the level of free iron in biological fluids. Serum transferrin level could be monitored for iron deficiency, hemochromatosis, and other iron overload disorders.

免疫原: A synthetic peptide near the N-terminus of human transferrin protein.

经过测试的应用: ELISA, WB, IHC, IP

推荐稀释比:

ELISA: 1:1000-1:5000

WB: 1:100-1:2000

IHC: 1:50-1:100

Concentration: 1 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG

纯化: Purified from ascites

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl and 50% glycerol


种属反应性: transferrin antibody recognizes transferrin of vertebrates

储存条件: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing.

Download Datasheet: 

蛋白质印迹:

Transferrin Protein



WB: Anti-Transferrin mAb

Purified His-tagged transferrin protein (货号:10192) was used for Western blot with anti transferrin monoclonal antibody (货号: 26050).