

TRANSFERRIN PAB

产品名称: Transferrin 兔多抗

货号: 21116

基因符号: TF

描述: Transferrin 兔多抗

背景: Transferrins are iron-binding blood plasma glycoproteins that control the level of free iron in biological fluids. Serum transferrin level could be monitored for iron deficiency, hemochromatosis, and other iron overload disorders.

免疫原: A synthetic peptide near the N-terminus of human transferrin protein.

经过测试的应用: ELISA, WB, IHC, IP

推荐稀释比:

ELISA: 1:1000-1:5000

WB: 1:100-1:2000

IHC: 1:100-1:500

Concentration: 0.2 mg/ml

种属反应性: Rabbit

形式: Liquid

克隆性: Polyclonal

亚型: IgG


纯化: Purified from serum

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol
Species

种属反应性: transferrin antibody recognizes transferrin of vertebrates

储存条件: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

蛋白质印迹:



WB: anti-transferrin pAb

NIH3T3 cell lysate was used for Western blot with anti-transferrin 兔多抗 (Cat # 21116).