

## CXCR4 MAB (ULOCUPLUMAB)

货号: 28101

产品全名: CXCR4(Ulocuplumab) 单抗

基因符

号CD184;D2S201E;FB22;HM89;HSY3RR;LAP-3;LAP3;LCR1;LESTR;NPY3R;NPYR;NPYRL;NPYY3R;WHIM;WHIMS

描述: CXCR4(ulocuplumab) 单抗

背景: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

经过测试的应用: Flow Cyt

推荐稀释比: Flow Cyt 1:100

种属反应性: Homo sapiens

亚型: IgG4

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: 人 CXCR4

成分: Lyophilized from sterile PBS, pH 7.4. 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at  $-80^{\circ}\text{C}$  (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 15  $\mu\text{g}/\text{mL}$  Anti-CXCR4 (ulocuplumab biosimilar) mAb (28101) on Expi293 cells transfected with Human CXCR4 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).