

## HUMAN UTS2R FULL LENGTH PROTEIN

货号: 11117

产品全名: 人 UTS2R Full Length 蛋白

规格: 10/50/100 µg

基因符号 GPR14; UR-2-R; UTR; UTR2

目标蛋白: UTS2R

**UNIPROT ID:** Q9UKP6

描述: 人 UTS2R full length protein-synthetic nanodisc

背景: High affinity receptor for urotensin-2 and urotensin-2B. The activity of this receptor is mediated by a G-protein that activate a phosphatidylinositol-calcium second messenger system.[UniProtKB/Swiss-Prot Function]

物种/宿主: HEK293

分子量: The human full length UTS2R protein has a MW of 41.9 kDa

**Formulation & Reconstitution:** Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.